

Experience new levels of control



GORE® TAG® Conformable Thoracic Stent Graft with ACTIVE CONTROL System



Combining the precise control of an intuitive delivery system with the exceptional conformability of the Conformable GORE® TAG® Thoracic Endoprosthesis, clinicians can achieve predictable patient outcomes.

- Staged deployment enables refinement of placement and angulation
- Optional angulation control at intermediate and full diameter
- Hemodynamic stability throughout the deployment
- Reduced profiles on 10 device sizes, with no changes to the stent graft itself



GORE® TAG® Conformable Thoracic Stent Graft with ACTIVE CONTROL System

Catalogue number	Intended aortic diameter (mm)	Proximal diameter (mm)	Distal diameter (mm)	Endoprosthesis length (cm)	Recommended GORE® Dryseal Flex Introducer Sheath size (fr)	Oversizing range (%)	Partially uncovered stent length (mm)
TGM212110	16–19.5	21	21	10	18	8–31	3
TGM262110	19.5–24 / 16–19.5	26	21	10	20	8–33	4
TGM262610	19.5–24	26	26	10	20	8–33	4
TGM282810	22–26	28	28	10	20	8–27	4
TGM282815	22–26	28	28	15	20	8–27	4
TGMR312610	24–29 / 19.5–24	31	26	10	20	7–33	4
TGMR313110	24–29	31	31	10	20	7–29	4
TGMR313115	24–29	31	31	15	20	7–29	4
TGMR313120	24–29	31	31	20	20	7–29	4
TGM343410	27–32	34	34	10	22	6–26	5
TGM343415	27–32	34	34	15	22	6–26	5
TGM343420	27–32	34	34	20	22	6–26	5
TGMR373710	29–34	37	37	10	22	9–28	5
TGMR373715	29–34	37	37	15	22	9–28	5
TGMR373720	29–34	37	37	20	22	9–28	5
TGMR404010	31–37	40	40	10	22	8–29	6
TGMR404015	31–37	40	40	15	22	8–29	6
TGMR404020	31–37	40	40	20	22	8–29	6
TGM454510	34–42	45	45	10	24	7–32	6.5
TGM454515	34–42	45	45	15	24	7–32	6.5
TGM454520	34–42	45	45	20	24	7–32	6.5

Learn more at goremedical.com/active-control



W. L. GORE & ASSOCIATES, INC.
Flagstaff, AZ 86004

+65.67332882 (Asia Pacific)
00800.6334.4673 (Europe)
800.437.8181 (United States)
928.779.2771 (United States)

goremedical.com

Products listed may not be available in all markets.

GORE, ACTIVE CONTROL, TAG and designs are trademarks of W. L. Gore & Associates.
© 2019 W. L. Gore & Associates, Inc. AY1110-EN1 AUGUST 2019

INDICATIONS FOR USE IN THE U.S.: The GORE® TAG® Conformable Thoracic Stent Graft is intended for endovascular repair of all lesions of the descending thoracic aorta, including: Isolated lesions in patients who have appropriate anatomy, including: adequate iliac / femoral access, aortic inner diameter in the range of 16–42 mm, ≥ 20 mm non-aneurysmal aorta proximal and distal to the lesion; Type B dissections in patients who have appropriate anatomy, including: adequate iliac / femoral access, ≥ 20 mm landing zone proximal to the primary entry tear; proximal extent of the landing zone must not be dissected, diameter at proximal extent of proximal landing zone in the range of 16–42 mm. **CONTRAINDICATIONS:** Patients with known sensitivities or allergies to the device materials; patients who have a condition that threatens to infect the graft. Refer to *Instructions for Use* at goremedical.com for a complete description of all warnings, precautions, and adverse events. **INDICATIONS FOR USE UNDER CE MARK:** The GORE® TAG® Conformable Thoracic Stent Graft is indicated for endovascular repair of the descending thoracic aorta. **CONTRAINDICATIONS:** Patients with known sensitivities or allergies to the device materials; patients with a systemic infection who may be at increased risk of endovascular graft infection. Refer to *Instructions for Use* at goremedical.com for a complete description of all warnings, precautions, and adverse events. Rx Only